

Connecting



from

heart

to

heart



SonoTrax Series Ultrasonic Pocket Doppler's superior sensitivity, high-fidelity sound, recording option and optional rechargeable battery ensure the doppler family to meet the routine examination requirements of obstetricians and midwives.

SonoTrax

Ultrasonic Pocket Doppler

Features

- Easy-to-use
- Compact design
- Interchangeable waterproof probe
- Built-in audio recorder
- Large backlight LCD display
- Accurate fetal heart detection with high-fidelity sound
- Support 2MHz, 3MHz, 4MHz, 8MHz waterproof probes
- Automatic and manual counting of FHR
- Rechargeable battery optional
- Automatic power-off

SonoTrax

Ultrasonic Pocket Doppler



Compact Li-ion Battery Charger



Backlight LCD



Built-in Recorder



Interchangeable Waterproof Probe



Specifications

Display: 45x25mm LCD display

Heart Rate Range: 50~210bpm

Resolution: 1bpm

Recording Time: 240s

Audio Output Power: 0.5W

Auto Power-off: 1 minute with no signal or operation

Size: 32 (Depth) x 85 (Width) x 138 (Height) mm

Weight: 290±5g (including one battery)

Model	Interchangeable Probes	FHR Display	Built-in Audio Recorder	Backlight	Waterproof Probe	Probe Auto-detection Display	Battery Status Display	Low Battery Indication	AC Recharge
<i>SonoTrax Lite</i>	Yes	No	No	No	Yes	No	No	Yes	No
<i>SonoTrax Basic</i>	Yes	Yes	No	No	Yes	Yes	No	Yes	No
<i>SonoTrax Basic A</i>	Yes	Yes	No	Yes	Yes	Yes	Yes	Yes	No
<i>SonoTrax Pro</i>	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	No
<i>SonoTrax II</i>	Yes	Yes	No	Yes	Yes	Yes	Yes	Yes	Yes
<i>SonoTrax II Pro</i>	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes	Yes
Optional accessories	Rechargeable Li-ion battery, charger, earphone, carry bag, 9V alkaline battery, adapter, Ni-MH battery, probe (2MHz, 3MHz, 4MHz, 8MHz)								



Edan Instruments, Inc.

3F-B, Nanshan Medical Equipments Park, Nanshan Rd 1019#, Shekou, Nanshan, Shenzhen, 518067 P.R. China
Tel: 86-755-2689 8326 Fax: 86-755-2689 8330 www.edan.com.cn Email: info@edan.com.cn



All rights reserved.

Features and specifications are subject to change without notice

ENG-CTG-SonoTrax-V1.4-20100310