



Excellence in Human Service  
And Technology

# MINIDOP ES-100VX

POCKET DOPPLER WITH GREAT SOUNDS

## FEATURES

- **Pocketable, Hand-held Size**
- **Great sensitivity**
- **Wide probe selection for user flexibility**  
Clinics: 2, 4, 5, 8 & 10MHz  
Surgical: 8, 10 & 20MHz
- **Convenient probe activation button**  
turn the unit ON and OFF

## CLINICAL APPLICATIONS

- **Blood flow detection**



## PROBE SELECTION

### — Standard type —



2MHz Probe



4, 5, 8 & 10 MHz Probe

### — Sterilizable type \*Amplifier required —



Small pencil type



Long pencil type

..and more

## SPECIFICATIONS

**Power Supply:** DC 9V, 9V alkaline square battery  
**Battery life:** Approx. 3 hours  
**Max speaker output:** 300 mW or more  
**Multi frequency:** 2, 4, 5, 8, 10 & 20MHz  
**Automatic shut-off:** No signal: 3min.  
**External output:** Headphone Connector: Speaker cut off

**Dimensions:** 78 (W) x 141 (L) x 27(H) mm  
(Probe holder not included)  
**Weight:** 240g (Battery and probe are included)  
**Accessories:** Carrying bag, Gel, Battery  
**Probe button:** Power ON/OFF  
**Electrical safety:** Conform to IEC60601-1

*\*Specifications to subject to change.*

### Hadeco, Inc.

2-7-11, Arima, Miyamae-ku, Kawasaki,  
Kanagawa, 216-0003, Japan  
TEL: +81 44 877 4361  
FAX: +81 44 855 7301  
<http://www.hadeco.co.jp/>



080-00125-07  
Printed in Japan